

RNA extraction

Reagents and materials

- Trizol
- Chloroform
- Eppendorf tubes
- Centrifuge
- Isopropanol
- 75% ethanol (etOH)
- Heating block
- Nanodrop

Trizol is highly toxic. Therefore, always handle it in the hood and with gloves and a lab coat, and use special waste containers.

Step 0

Pre-cool down the centrifuge (4°C) and turn on the heating block (60°C).

Step 1

Aspirate the media of the cells and add 400µL Trizol per well (from the 6cm plate).

Step 2

Pipette up and down and transfer all to an Eppendorf tube. Wait 5 minutes at RT.

Step 3

Add 80µL of chloroform. Invert the tubes five times and incubate 2-3 minutes at RT.

Step 4

Centrifuge 15 minutes at 12.000g (4°C). Transfer the aqueous phase (upper part) into a fresh Eppendorf tube. BE CAREFUL TO NOT TAKE ANYTHING FROM OTHER PHASES!

Step 5

Add 200µL of isopropanol and incubate 10 minutes at 4°C.

Step 6

Centrifuge 10 minutes at 12.000g (4°C). Discard the supernatant with a pipette and resuspend the pellet with 400µL of 75% etOH.

Step 7

Mix well and centrifuge 5 minutes at 7.500g (4°C). Discard the supernatant with a pipette and air dry the pellet for 10 minutes at 60°C in the heating block.

Step 8

Resuspend the pellet in 20-50µL of RNase free water, depending on the size of the pellet.

Step 9

Quantify the RNA using the Nanodrop and store at -80°C until cDNA synthesis.

cDNA synthesisReagents and materials

- SCRIPT cDNA synthesis kit (#PCR-511L) from Jena Bioscience.
- Eppendorf tubes
- Heating block

Step 1

Add the following components to a nuclease-free microtube. Pipette the mix on ice and mix gently by pipetting up and down. Use the document “cDNA synthesis calculations template”.

Component	Stock concentration	Final concentration	10μL assay
RNAse-free water	-	-	Fill up to 10 μ L
RNA template	-	1 μ g	1 μ g (calculate depending on RNA concentration)
Primer (random primer)	100 μ M	50pmol (100ng)	0.25 μ L
SCRIPT RT Buffer complete	5x	1x	2 μ L
dNTP Mix	10mM each	500 μ M each	0.5 μ L
DTT stock solution	100mM	5mM	0.5 μ L
SCRIPT Reverse Transcriptase	200 units/ μ L	100 units	0.25 μ L

Step 2

First strand cDNA synthesis:

- Incubate at 42°C for 10 minutes.
- Incubate at 50°C for 45 minutes.

Step 3

Inactivation of the Reverse Transcriptase:

- Heat the mixture to 70°C for 10 minutes.

Step 4

Store at -20°C in a 1:20 dilution (10 μ L cDNA + 190 μ L H₂O).

qPCR

Reagents and materials

- SsoAdvanced™ Universal SYBR® Green Supermix (BioRad, #1725275)
- MicroAmp Fast Optical 48-Well Reaction Plate
- MicroAmp™ 48-Well Optical Adhesive Film (25 films)
- Primers
- cDNA
- H₂O
- Centrifuge
- StepOne qPCR machine

Step 1

Mix to a final volume of 5µL/well:

- 2,5µL Power SYBR™ Green PCR Master Mix (2x)
- 0,15µL forward primer (10µM)
- 0,15µL reverse primer (10µM)
- 2µL cDNA
- 0,2 µL H₂O

Use the document “qPCR calculations and plate design template”.

Step 2

Seal the plate with the adhesive film.

Step 3

Spin down the plate in the centrifuge and start the following PCR program.

Holding stage:

- 10 min at 95°C

Cycling Stage (40 cycles):

- 15 sec 95°C
- 1 min 60°C

Melt curve stage:

- 15 sec 95°C
- 1 min 60°C
- 60°C → 95°C

